



zwitter ion of glycine:




Net charge = 0.




zwitter ion (at pI)

(pH > pI) charge = -1




zwitter ion (at pI)

(pH < pI) charge = -1




low pH



zwitter ion

$\xrightleftharpoons{K_{a2}}$



high pH

Aim: Determination of pI of Glycine using pH-metry.

Apparatus Required: pH-meter, 100 ml beaker, burette, pipette, flasks, glass rods.

Chemicals required: NaOH solution (0.1 M), glycine, distilled water, dil. HCl.

Principle: Potentiometry deals with the measurement of potential differences between two electrodes, which are combined to form an electrochemical cell. Glycine is the simplest amino acid. In aqueous solution, glycine is amphoteric; at low pH it can be protonated and at high pH it loses a proton. Zwitter ions have equal number of positive and negative charges and are electrically neutral overall.

Isoelectric point: Isoelectric point (pI) is the pH at which the amino acid has an overall zero charge and exist in the zwitter ionic form. The pI can also be understood as the pH at which the amino acid doesn't migrate in an electric field.

The pI is given by the average of the two  $pK_a$ . where  $K_a$  is the dissociation constant for loss of proton.

• In the solutions more acidic than pI;

- (1) the  $\text{COO}^-$  in the zwitter ion accepts a proton.
- (2) the amino acid has positive charge.

• In the solutions more basic than pI,

- (1) the  $\text{NH}_3^+$  in the zwitter ion loses a proton.
- (2) the amino acid has a negative charge.



Observations :

Concentration of NaOH = 0.1N.

Volume is in (ml)

Vol NaOH	pH	Vol NaOH	pH	Vol NaOH	pH	Vol NaOH	pH
0	1.86	20	2.66	31.1	4.07	46.5	9.53
0.5	1.86	20.5	2.69	31.3	4.21	47	9.56
1	1.87	21	2.73	31.4	4.35	47.5	9.57
1.5	1.89	21.5	2.75	31.5	4.41	48	9.59
2	1.91	22	2.78	31.6	4.68	48.5	9.62
2.5	1.93	22.5	2.81	31.7	5.75	49	9.65
3	1.95	23	2.85	31.9	6.98	49.5	9.7
3.5	1.96	23.5	2.89	32	7.21	50	9.72
4	1.98	24	2.9	32.1	7.52	50.5	9.74
4.5	2	24.5	2.93	32.2	7.71	51	9.75
5	2.01	25	2.97	32.3	7.8	51.5	9.79
5.5	2.03	25.5	3.02	32.4	7.89	52	9.83
6	2.04	25.6	3.01	32.5	7.95	52.5	9.86
6.5	2.06	25.7	3.02	32.6	8.04	53	9.89
7	2.08	25.8	3.04	33	8.1	53.5	9.91
7.5	2.1	25.9	3.05	33.5	8.2	54	9.93
8	2.11	26	3.06	34	8.36	54.5	9.97
8.5	2.13	26.5	3.1	34.5	8.49	55	10
9	2.15	27	3.15	35	8.59	55.5	10.04
9.5	2.16	27.5	3.21	35.6	8.66	56	10.07
10	2.19	27.6	3.22	36	8.77	56.5	10.1
10.5	2.21	27.7	3.24	36.5	8.79	57	10.13
11	2.23	27.8	3.24	37	8.86	57.5	10.16
11.5	2.25	27.9	3.25	37.5	8.91	58	10.19
12	2.27	28	3.27	38	8.96	58.5	10.23
12.5	2.3	28.5	3.34	38.5	9	59	10.25

13	2.33	28.6	3.36	39	9.04	59.5	10.33
13.5	2.34	28.8	3.38	39.5	9.08	60	10.36
14	2.37	28.9	3.4	40	9.11	60.5	10.4
14.5	2.4	29	3.51	40.5	9.15	61	10.44
15	2.43	29.5	3.53	41	9.2	61.5	10.5
15.5	2.45	29.6	3.56	41.5	9.23	62	10.55
16	2.5	29.7	3.58	42	9.27	62.5	10.6
16.5	2.51	29.8	3.6	42.5	9.29	63	10.67
17	2.53	29.9	3.64	43	9.32	63.5	10.71
17.5	2.55	30	3.67	43.5	9.36	64	10.76
18	2.57	30.1	3.78	44	9.38	64.5	10.83
18.5	2.6	30.5	3.82	44.5	9.4	65	10.87
19	2.63	30.9	3.93	45	9.44	65.5	10.94
19.5	2.65	31	3.99	45.5	9.45	66	10.99



Procedure :

- 1) Take 50ml of  $\text{NaOH}(0.1\text{M})$  solution in a 50ml burette and adjust zero reading.
- 2) Pipette out 25ml of the given solution of amino acid in a 250ml plastic beaker and add 25ml of distilled water to the amino acid solution using a pipette.
- 3) Insert the cleaned pH electrode into the beaker solution and record the initial pH of the solution.
- 4) Don't remove the electrode from the beaker till the end of experiment.
- 5) Add  $\text{NaOH}$  in 0.5ml increments from the burette and stir the solution and mix it well.
- 6) Record the corresponding pH values until the pH starts increasing drastically. At this time, add 0.1 ml increments of  $\text{NaOH}$  till the pH stabilises around 8.
- 7) After reaching pH 8, continue adding 0.5ml increments of  $\text{NaOH}$  till you reach pH 11.
- 8) Plot the graph of pH vs volume of  $\text{NaOH}$  solution.
- 9) The two almost horizontal parts of the graph give the values of  $\text{pK}_1$  and  $\text{pK}_2$  of glycine. Use mid points of these regions to get the values.
- 10) The average of these values ( $\text{pK}_1$  and  $\text{pK}_2$ ) gives the pI of glycine.

Teacher's Signature \_\_\_\_\_

from the graph:

$$pK_{a1} = 2.45$$

$$pK_{a2} = 9.65$$

$$pI = \frac{1}{2}(pK_{a1} + pK_{a2}) = \frac{1}{2}(2.45 + 9.65) = \frac{12.10}{2} = \underline{\underline{6.05}}$$

Result:

- 1) pI of glycine was determined using pH-metry.
- 2) pI of glycine = 6.05.

Precautions:

- 1) Handle the glassware carefully.
- 2) pH meter should be properly calibrated.
- 3) Sensitivity of electrode should be taken care of.



Scale: 1 cm = 5 ml on X axis  
1 cm = 1 unit pH on Y axis

